



DICOM Conformance Statement for WEB DICOM Viewer

MedDream

Version 2.01

2009-10-08

MedDream DICOM v3.0 conformance summary

Introduction

This section is an abbreviated DICOM conformance statement for MedDream DICOM viewer.

Supported transfer syntaxes (Reading)

The Transfer Syntax UID is in the file DICOM Tag field (0002,0010).

| Uncompressed Transfer Syntax | Description |
|------------------------------|--|
| 1.2.840.10008.1.2 | Implicit VR - Little Endian |
| 1.2.840.10008.1.2.1 | Explicit VR - Little Endian |
| 1.2.840.10008.1.2.1.99 | Deflated Explicit VR - Little Endian |
| 1.2.840.10008.1.2.2 | Explicit VR - Big Endian |
| 1.2.840.113619.5.2 | Implicit VR - Big Endian (G.E Private) |
| RLE Transfer Syntax | |
| 1.2.840.10008.1.2.5 | Run Length Encoding, Lossless |
| JPEG Transfer Syntax | |
| 1.2.840.10008.1.2.4.50 | JPEG Baseline (Process 1) |
| 1.2.840.10008.1.2.4.51 | JPEG Extended (Process 2 & 4) |
| 1.2.840.10008.1.2.4.57 | JPEG Lossless, Non-Hierarchical (Process 14) |
| 1.2.840.10008.1.2.4.70 | JPEG Lossless, Hierarchical, First-Order Prediction (Process 14) |
| 1.2.840.10008.1.2.4.90 | JPEG 2000 Image Compression (Lossless Only) |
| 1.2.840.10008.1.2.4.91 | JPEG 2000 Image Compression |
| 1.2.840.10008.1.2.4.100 | MPEG2 Main Profile @ Main Level |

Supported "Photometric Interpretation" pixel format (Reading)

The Photometric Interpretation UID is in the file DICOM Tag field (0028,0004).

| Photometric Interpretation pixel format | Description |
|---|--|
| MONOCHROME1 | grey level image description (high values=dark, low values=bright) |
| MONOCHROME2 | grey level image description (high values=bright, low values=dark) |
| PALETTE COLOR | pseudo color image description |
| RGB | true color image description |
| YBR_FULL | true color image description |
| YBR_FULL_422 | true color image description |

Supported 'Bits Allocated' values (Reading)

The Bits Allocated value is in the file DICOM Tag field (0020,0100).

| Classical values | Description |
|------------------|-------------|
|------------------|-------------|

8, 12, 16 12 means that 4 pixels are stored in 3 'short int'

| Unusual values | Description |
|----------------|-------------|
|----------------|-------------|

24 Some ACR-NEMA RGB files came with 'Bits Allocated' = 24 and "Samples Per Pixel" = 1, or with no "Samples Per Pixel"

32 Some ACR-NEMA files, "Bits Allocated" = 32