

Quick Install Guide for MedDream DICOM Viewer 7.5.2: Standard DB on PostgreSQL integration

Important: please note that instructions above are for Windows only.

1. Install PostgreSQL if it is not installed:

<https://get.enterprisedb.com/postgresql/postgresql-10.8-4-windows-x64.exe>

During the PostgreSQL server configuration remember the port and superuser password;

2. Install MedDream:

<https://www.softneta.com/files/meddreamviewer/752/MDPACK-7.5.2.exe>

3. Install MedDream DB if it is not installed yet for previous versions:

<https://www.softneta.com/files/meddreamviewer/752/MedDreamViewer-7.5.2-DB-PostgreSQL.zip>

3.1. Unzip MedDreamViewer-7.5.2-DB-PostgreSQL.zip to C:\MDPACK\MedDreamViewer-7.5.2-DB-PostgreSQL

3.2. Install MedDream DB to MySQL by executing C:\MDPACK\MedDreamViewer-7.5.2-DB-PostgreSQL\PostgreSQL-DB\meddream-db-create.bat

You will be asked to enter a new password for the user "meddream". Remember it.

Immediately after that you will be asked to enter the current password for the superuser.

3.3. Insert test data to MedDream DB by executing C:\MDPACK\MedDreamViewer-7.5.2-DB-PostgreSQL\PostgreSQL-DB\meddream-db-insert.bat

you will be asked to enter the password for the superuser.

4. Configure MedDream:

4.1. Copy and replace c:\MDPACK\MedDreamViewer-7.5.2-DB-PostgreSQL\MDPACK folder to C:\MDPACK

4.2. Configure meddream DB by entering user "meddream" password in C:\MDPACK\MedDream\application.properties file:

com.softneta.meddream.pacs.configurations[0].password={Password} #line 418

5. Restart windows services (Run in CMD line: services.msc):

5.1. Restart MedDream

You can connect to MedDream with Logins:

http://127.0.0.1/

User name: demo

Pass: demo

6. URL Integration:

<http://127.0.0.1/?study=1.2.826.0.1.3680043.8.1055.1.20160224130121426.315299872.5643427>