

Function	Description	v8.2.0
Specific VET Functions		
Cobb angle	Angle measurement between two lines.	+
Tibial Plateau Angle (TPA)	Measure the angle between the tibial plateau slope line and the line that is perpendicular to the tibia axis line.	+
Norberg Angle	Angle measurement to evaluate canine hips.	+
Verbal Heart Scale (VHS)	Measure heart size and provide an accurate assessment of true cardiac enlargement.	+
Hip distraction index (Hip DI)	Measurement allows evaluation of the laxity of the hip joint.	+
Tibial tuberosity advancement (TTA)	Measurement allows for estimating the cage width of tibial tuberosity advancement (TTA) in the surgical stabilization of the dog's knee joint.	+
<p>This specification file only describes only Specific MedDream DICOM VET Viewer veterinary functions. Please read the Specification document for a full list of MedDream DICOM Viewer specifications.</p>		